



AURA™



Polycentric Knee

The Aura polycentric knee was designed to offer stability and relief to low-impact users. Featuring a unique push-button lock and simple 4mm adjustments, the waterproof Aura goes beyond mere physical support. It can open doors for mobility goals and newfound peace of mind.



FEATURES:



REFINED ADJUSTMENTS

Like all College Park knees, Aura adjustments are quick and easy. Clinicians can adjust stance flexion, extension assist, and more with a 4mm hex key.



UNIQUE LOCK AND GEOMETRIC STABILITY

The Aura's proximal center of rotation was designed to offer inherent stability. Users can easily lock and unlock through the simple push of a button. This unique method may improve the user experience and enable usage in water.



WATERPROOF

Designed and tested for corrosion resistance, the Aura has been approved for submersion in the ocean and chlorinated pools. It was built to perform well and maintain high quality in different wet environments.

TECHNICAL SPECIFICATIONS

WEIGHT	OVERALL HEIGHT	WARRANTY	WATER CAPABILITY
665 g	161 mm	2 years	 Waterproof
FLEXION RANGE	DISTAL ATTACHMENT	WEIGHT LIMIT	DOME TO KNEE CENTER
135°	30 mm Tube Clamp	136 kg	17 mm





AURA PART NUMBER

APKA 150

